

Conductivity Meter
EC-3210

Conductivity Meter
EC-3210



Features:

- Microprocessor-Based
- Selectable Cell Constants
- Selectable Temperature Coefficients
- Auto Ranging Capability
- Automatic or Manual Temperature Compensation
- Large LCD Display

Specifications:

Range Conductivity	0.05~ 2x10 ⁵ μ S/cm
Temp	0~ 100°C
Accuracy Conductivity	±1.0% Full Scale
Temp	±1°C
Temp Compensation	0~ 40°C, Automatic or Manual
Temp Coefficient	0~ 5%/°C
Dimension	210x205x65mm
Power	DC9V or 220VAC adapter, 50Hz
Weight	1.5kg

PHYLO FOR TECHNOLOGY

P.iva 03176751208
VIA EDOARDO FERRAVILLA,
19
40127-BOLOGNA

ITALY

Tel : + 390 51 4126484
Fax: + 390 51 9913348