

Cryo Tk - 50

50L

PHYLO FOR TECHNOLOGY

P.iva 03176751208

VIA EDOARDO FERRAVILLA, 19

40127-BOLOGNA

ITALY

Tel :+ 390 51 4126484

Fax:+ 390 51 9913348

Website: www.phylo-tec.com
E-Mail: info@phylo-tec.com

Cryogenic Tanks



Cryo Tk-50

50L

Characteristic:

The product is made of stainless steel and equipped with automatic pressure-increasing system, available to form pressure to transfuse liquid continuously by itself. With the characteristics of high efficiency and continuance, it can mainly be used in transferring transportation. It may also safely and conveniently be applied to liquid supplement to biological Dewars of liquid nitrogen, liquid oxygen and liquid argon, or used together with electronic equipment.

Technical parameters:

- Geometric Volume: 55L
- Available Volume: 50L
- Dimensions (mm): 450LX450WX1000H
- Empty Weight: 55Kg
- Full Weight: 95.4Kg
- Working Pressure: ≤ 0.2 Mpa
- Daily Evaporating Rate: < 0.02 L/day
- Transfusion Amount: ≥ 4 L/Min